



MASSAGE GUN

LS-016



Let's start your experience



*The illustrations of products, accessories, user interfaces, etc. in this manual are schematic diagrams for reference only. Due to product updates and upgrades, the actual product may slightly differ from the schematic diagrams. Please refer to the actual product.

Operation Instructions

Power On/Off: Power On: Press and hold the power button for 1.5 seconds to turn on the device
Gear Shift: The device defaults to "Gear 1" when powered on. Short press the power button to switch gears (initial gear is Gear 1, LED1 white light display). Short press to enter Gear 2 (LED2 white light display), Short press to enter Gear 3 (LED3 white light display).
After Gear 3, short press to return to Gear 1, and so on in a cycle.
Blockage Protection & Restart: If excessive pressure is applied for more than 3 seconds, blockage protection will be triggered. The device will pause for 3 seconds and then attempt to resume normal speed (motor vibration may be felt). If it fails to recover after 5 attempts, it will shut down.
Timed Shutdown: The device will automatically shut down after 15 minutes of default operation. Any operation within 15 minutes will not affect the running time.

- Select a suitable massage head, switch gears according to your personal experience, and it is recommended to start from a low gear for first-time use
- Gently press the massage head against the muscle area that needs relaxation, slide it back and forth slowly. If there is a local area with extra soreness, you can apply appropriate pressure on the tender point and stay for 15-30 seconds
- The fascia gun itself has impact force, so the pressure applied should not be too strong to avoid injury to muscles and bones
- It is recommended to use it on areas with plump muscles, and do not stay on the same area for more than 60 seconds
- The recommended usage duration is 30 minutes. If you feel soreness, you can adjust the massage time appropriately, with the total usage duration not exceeding 60 minutes

Battery Display: In standby mode, press the power button to display battery level (when battery is ≥ 100%, LED1, LED2, LED3 white lights are on; when battery is <70%, LED1, LED2 white lights are on; when battery is <30%, only LED1 white light is on; when battery is <10%, LED1 red light is on)
In power-on mode, when battery is <10%, LED1, LED2, LED3 white lights flash simultaneously, the device will automatically shut down when battery voltage is below 6V


Note: LED1 is a dual-color light (red/white)
Charging: When battery is <30%, LED1 red light flashes; when battery is >30%, LED1 red light turns off and white light flashes; when battery is >50%, LED1 white light stays on and LED2 white light flashes; when battery is >70%, LED1, LED2 white lights stay on and LED3 white light flashes; when fully charged, LED1, LED2, LED3 stay on
Note: The device cannot be started while charging, and buttons are ineffective

Four Massage Heads

- Original Metal Flat Head: Penetrates deeply to relax muscles with strong force feedback, suitable for thick muscle groups such as the legs
- Silicone Shock-absorbing Head: Delivers gentle impact, suitable for muscle groups such as the pectoralis major and latissimus dorsi
- Silicone Cone Head: Targets and relieves tender points, ideal for impacting deep tissues such as the palms and soles of the feet
- Silicone U-shaped Head: Designed for muscle soreness areas and special positions such as the muscles on both sides of the spine

⚠ WARNING

- When holding the device, please keep your fingers away from the percussion head to avoid pinching injury
- Install and remove the massage head only when the device is turned off. Ensure the massage head is securely installed before starting the device
- Do not adjust or remove the massage head during use to avoid pinching injury



Watch your fingers, beware of pinching injury

Product Basic Parameters

Product Name	Fascia Gun		
Model	LS-016	Rated Charging Voltage	7.4V ~
Battery Capacity	600mAh	Rated Input Power	10W

Daily Maintenance



Avoid direct sunlight



Avoid high temperature and humidity



Avoid proximity to fire and smoke



Avoid dusty areas

- Keep away from hard objects to avoid collision and scratch
- For external cleaning of the device, use a wet wipe or soft cloth
- If not used for a long time, it is recommended to charge once a month and keep the battery fully charged when in use

Troubleshooting

If the product does not work properly, please check for the following conditions; if the problem still cannot be solved, please return the product to the after-sales service for repair

Condition	Possible Causes	Solutions
Cannot start	Low battery	Charge the device
	Device not turned on	Press and hold the power button to turn on
	Device still connected to household power socket	Remove the charging cable from the device
Cannot charge or no charging light display	Charging cable not properly inserted into the power socket	Replug the charging cable
	Incorrect charging cable used	Use a Type-C data cable

With you always
"Wish you a relaxed and happy mood every day"



Product Disposal Notes

Please support the recycling of rechargeable batteries (lithium-ion batteries)

	Danger The battery used in this product is a dedicated rechargeable battery and must not be used in other products. Do not charge the rechargeable battery after it has been removed. • Do not throw it into fire or expose it to excessive heat • Do not pierce it with nails, subject it to heavy impact, or disassemble and modify it • Do not directly connect the positive and negative terminals of the battery with metal • Do not store it together with necklaces, hairpins, etc. • Do not charge, use, or place it in high-temperature places such as near stoves or in direct sunlight All of the above situations can cause heating, fire, explosion, or rupture.
	Warning Keep the removed rechargeable battery out of the reach of infants and young children to prevent accidental swallowing • If accidental swallowing occurs, seek medical attention immediately • If the battery leaks, do not touch the liquid with bare hands and handle it as follows • When battery fluid gets into the eyes, it may cause blindness. Do not rub the eyes directly with hands. Rinse immediately with clean water and seek medical attention. If battery fluid gets on skin or clothing, it may cause skin inflammation or injury. Rinse with clean water and seek medical attention if discomfort persists • Disassembly is only allowed when the product is being discarded. Failure to comply may result in fire, electric shock, injury, or other accidents. For repairs (except battery replacement), please consult the store where you purchased the product

The lithium-ion battery used in this product is a valuable recyclable resource. Please hand it over to a professional institution that handles rechargeable battery recycling and regeneration when the product is discarded after use.

- Before discarding the device, the battery must be removed from it
- When removing the battery, the device must be powered off
- The battery must be disposed of safely

Removal Steps:

1. Use a tool to open the outer casing
2. Unscrew the screws on the circuit board
3. Unplug the battery terminal connector
4. Remove the battery

Precautions

- Please read the instruction manual carefully before using the product and keep it properly
- Use and charge the product in an environment of approximately 0-40 °C
- Keep away from water, fire, and corrosive environments to avoid affecting battery performance
- Do not repair, disassemble, or modify the product by yourself
- Do not use the product while it is charging
- After charging is complete, please unplug the charger from the product in a timely manner

Usage Prohibitions

- Do not use this product on skeletal areas with less muscle coverage, such as the head, elbows, knees, joints, and spine
- Do not use this product on the lateral or anterior side of the neck to avoid injury to the cervical arteries, which may cause fainting or cardiac arrest
- Do not use this product on areas with dense and fragile blood vessels, such as the armpits, inner upper arms, inner clavicle, inner elbow joints, popliteal fossa, abdominal kidney area, and genitals
- Patients with urinary system diseases such as kidney stones, ureteral stones, bladder stones, urethral stones, and urinary tract infections are prohibited from using this product below the abdomen
- Do not use this product on areas with bruising, body pain, surface wounds, local fever, redness and swelling, acute inflammation, muscle spasms, muscle tears, acute sprains, or contusions
- After strenuous exercise, if symptoms such as general fever, nausea, vomiting, or muscle pain and weakness occur, stop using this product immediately and seek medical attention
- Patients with the following diseases are prohibited from using this product: cardiovascular and cerebrovascular diseases, blood diseases, coagulation dysfunction, thrombocytopenic purpura, neck plaques, varicose veins, cancer, tumors, and mental illnesses
- Patients with fractures or internal steel nails are prohibited from using this product
- Do not use this product during precision operations or driving
- Do not use this product continuously for a long time to prevent muscle injury caused by excessive massage

If you have any of the following conditions, consult a doctor before use, otherwise accidents or physical discomfort may occur; if you meet any of the following conditions, consult a surgeon or attending physician in a Grade A tertiary hospital before use, otherwise accidents or physical discomfort may occur; minors and frail elderly people should use this product under the supervision and guidance of family guardians; if you have congenital diseases, do not use it privately to avoid accidents, and you can only use it with a doctor's permission and accompanied by others!

- Patients with heart or cranial nerve abnormalities
- Rhabdomyolysis syndrome
- Those with a cardiac pacemaker or stent implanted
- Those who have undergone surgical treatment within 12 months
- Those with abnormalities of the spine, joints, or skeleton, such as osteoporosis, cervical spondylosis, lumbar disc herniation, and other orthopedic diseases
- Those with abnormal limb sensation or mobility impairment caused by disease or age
- Those with complications caused by other diseases, such as psychiatric disorders and retinal damage

Warnings

- This product is a personal massage device and cannot replace drugs or medical devices
- When the device is running, do not place fingers, hair, or other body parts near the massage head and device interface to avoid injury
- Do not block or insert foreign objects into the massage head mounting slot or charging port, as this may cause burns, electric shock, or fire due to short circuit
- Do not immerse this product in water or other liquids
- Do not use this product in humid environments
- Do not plug or unplug the charger with wet hands, as this may cause electric shock or other accidents
- Do not disassemble this product by yourself
- Do not replace the built-in battery by yourself
- Please use in the normal operating environment: 0 °C~40 °C; relative humidity <80%RH, no condensation; atmospheric pressure: 50kpa-106kpa
- Please choose normal transportation and storage conditions: ambient temperature range: -20 °C to +55 °C; relative humidity range <95%RH, no condensation; atmospheric pressure: 50kpa-106kpa
- If you feel any discomfort during use, please stop using immediately

Prohibited Use in the Following Cases

- | | | |
|---|--|--|
| • Patients with malignant tumors | • Pregnant women and postpartum women | • Cardiac pacemaker users |
| • Near the human heart | • Patients with bleeding tendency | • Patients with acute suppurative inflammation |
| • Patients with metal implants in the body | • Patients with spastic paralysis | • Patients with coagulation dysfunction |
| • Those in the spasm stage after cerebrovascular accident | • Those with severe cognitive impairment | |
| • Those with abnormal nerve irritability | • Those with impaired consciousness and limb or skeletal joint deformities | |

Attention

- Do not use high speed and high pressure on the head or any skeletal parts of the body
- Frequent massage on the same area may cause injury; stop using immediately if pain or discomfort occurs
- Keep fingers, hair, or any body hair away from the moving parts of the device to prevent pinching or injury
- Do not misuse this product
- People with contraindicated symptoms should use this product only with a doctor's permission
- When starting the product, ensure there are no obstacles in front of the massage head. If it pauses due to external resistance, try restarting it
- This product is not waterproof; do not immerse it in water or other liquids
- Do not disassemble the entire product. If accidental disassembly causes malfunction, contact after-sales service promptly (all resulting costs will be borne by the consumer)
- Prolonged continuous operation may cause the motor surface temperature to rise; use it reasonably

Other Abnormalities and Malfunctions

- Avoid dropping, shaking, or striking the product. Strong impacts may damage components and cause malfunction. If any of these occur, contact after-sales service for handling promptly.
- If the device exhibits excessive temperature, smoke, buzzing, deformation, swelling, or battery leakage during use, stop using it immediately and contact after-sales service.

Recommendations

When using for the first time or after long-term storage, the battery capacity may be depleted; please charge the device according to the instruction manual. When discarding this product, please dispose of it in accordance with the relevant environmental protection regulations of your city, otherwise it may cause environmental pollution. Recommended for users aged 16 and above.

Note: Do not use or store the product in the following locations

Places exposed to direct sunlight	Places with water
Near flames or fire	Places with strong electromagnetic fields
High-temperature and humid places	Places with excessive dust