


# KNEE MASSAGER

## Take care of family health



Product instructions



Thanks very much for choosing this product  
Please read this manual carefully before use ,  
and keep it for future reference.

### Product Features

- 1.Built-in 20 red light therapy lamps/air pressure  
massage/hot compress/vibration.
- 2.The 4-sided three-dimensional wrapping design  
takes care of your knee joints.
- 3.The product is equipped with a 2000mAh lithium  
battery, which is convenient to use.

Product Specification  
Rated Voltage: DC 5V  
Rated Current: 1A  
Rated Power: 5W  
Product Net Weight: 0.8Kg  
Product Size: 10\*6\*4.5 Inches  
Auto Shutoff Time: 15/20/25/30Mins  
Heating Temperature: Low Level-115°F; Medium Level-  
125°F;High Level-130°F  
Executive Standards: GB4706.1-2005, GB4706.10-2008

### Product Functions and Configurations



### User's Guidance

- 1.In order to achieve a better effect, please try to tie the  
adjustable strap of the machine tightly onto your knee.
- 2.Long press the "On/Off" button, the machine starts to  
work. After the product is turned on, the default functions  
are: low-level temperature, low-level kneading pressure,  
low-level vibration, and red light on.
- 3.After the product is turned on, short press the "Heating"  
button to adjust the heating temperature (low-H1, medium  
-H2, high-H3 and heating off-HF), the display on screen will  
be updated at the corresponding temperature position with  
a beep prompt;
- 4.Short press the "Intensity" button to adjust airbag  
kneading strength (weak-L1, medium-L2, strong-L3), the  
display on screen will be updated at the corresponding  
intensity position with a beep prompt;
- 5.Short press the "Mode" button to adjust the massage  
modes (vibration/air pressure/heating-P1, vibration-P2, air  
pressure-P3), the display on screen will be updated at the  
corresponding mode position with a beep prompt;
- 6.Short press the "Red light" button to turn on the red light  
physiotherapy function, and the display on screen will be in  
the corresponding red light position with a beep sound.
- 7.Short press the "Time" button to adjust the auto shutoff  
time(15mins-20mins-25mins-30mins), and the display on  
screen will be updated in the corresponding time position  
with a beep sound.
- 8.The continuous using time better not exceed 30 minutes  
for each time.

### Product Charging

- 1.Please use a DC 5V adapter, insert the USB end of the charging  
cable that comes with the device into the adapter, then insert the  
other end of the cable into the power jack of the machine, the  
display will show the charging status;
- 1)There are 5 groups of blue battery display icons at the bottom of  
the display screen, if all 5 groups of icons keep flashing, it means  
the machine starts charging and the current power is less than 33%;
- 2)If the charging completes 33%, the outer four groups of icons will  
keep flashing while the middle one is still on;
- 3)If the charging completes 66%, the outer two groups of icons will  
keep flashing while the middle three are still on;
- 4)When the charging is complete, all five groups of icons will on.
- 2.When using the machine, if the outer 2 group of the battery icon  
flashes, please stop using it and charge it immediately.
- 3.When charging the machine, please use an adapter with DC 5V,  
1A or more than 1A. Of course, you can also use the computer USB  
interface to charge the machine.

### Environmental Protection and Resource Recovery

This product uses lithium batteries, please be sure to dispose of the  
old batteries at the designated place in accordance with the relevant  
local regulations.

### How to remove the built-in rechargeable battery

- This product contains a rechargeable battery, please be sure to remove  
the battery out from the product before consumers decide to dispose  
of the product.
- 1.Disassemble the product and be sure to disconnect the power supply  
before removing the battery.
- 2.Follow the instructions to remove the battery step by step.
- 3.The battery must be discarded in the local recycling and processing  
unit regulated by the government.

### Precautions before use

- 1.Never go to sleep during the use of this product (there is a risk of  
burns).
- 2.Patients with acute swelling or bleeding of the knee should consult  
a doctor before use.
- 3.Patients with rheumatoid arthritis should consult a doctor before use.
- 4.Patients with suspected edema, bacterial infection and other  
special inflammations caused by acute trauma should consult a doctor  
before use.
- 5.Pregnant women, infants and children, patients with skin intuition  
disorders or skin abnormalities are prohibited from using.
- 6.This product is designed as a household electronic product and is not  
suitable for diagnosing and treating diseases.
- 7.Please stop using the product immediately when the pain increases  
after using.
- 8.When using this product, if you feel that the kneading intensity is too  
strong, you can long press the "Intensity" button to release the air  
urgently.
- 9.In case of any fault, please ask professional and technical personnel  
for maintenance.
- 10.If the product is not in use for a long time, please put the product  
and the power adapter in the packaging box and store them in a dry  
and ventilated place. Do not dismantle and press for storage and keep  
them away from moisture and fire.
- 11.If it' s necessary to clean the product, please do not immerse the  
product in water or machine wash. Please use a soft brush or a clean  
cloth soaked with water or neutral detergent to clean it after the power  
is turned off. Only after the product is fully dry, it can be used with  
electricity.

### Safety Warning

Lithium battery warning  
This device contains a fixed lithium battery that must be recycled or  
disposed of properly.  
Contact your local waste authority for battery recycling and disposal.

### Notice

It needs to be charged in time to avoid excessive discharge when the  
power is too low.  
It should be charged every other month to avoid battery aging if it is  
not used for a long time.

### Common Problems and Solutions

#### Troubleshooting

If you are unable to solve the most common problems encountered  
when using this product with the information below, please consult  
the customer service center.  
If you have a problem not mentioned in the list, you should try the  
following steps first:  
1.Turn off the power of the product and unplug the power plug from  
the socket;  
2.Keep the product turned off for 60 minutes;  
3.Re-insert the plug into the power socket, and turn on the product  
after charging is complete.

| Problems                 | Possible Reasons                                                                                                 | Solutions                                                                                                                                      |
|--------------------------|------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Product does<br>not work | Charging product<br>does not work                                                                                | This is the safe protection<br>design. Unplug the charging<br>power after charging, then you<br>can press" on/off" to start for<br>normal use. |
|                          | The product just<br>completed 15 minutes<br>massage and<br>automatically switched<br>to standby mode             | One cycle is 15 mins,if you<br>need to massage again,please<br>press" Time" to select massage<br>time(15-20-25-30mins)                         |
|                          | Continuous use of this<br>product for 45mins<br>causes the thermal<br>protection device to<br>cut off the power. | Cool for 60<br>minutes before reuse.                                                                                                           |
| Probdunct is<br>heat     | Heating function on                                                                                              | Turn off the heating function.                                                                                                                 |
|                          | Use this product<br>continuously for 45<br>minutes                                                               | Turn off the product and<br>cool for 60 minutes.                                                                                               |
| Product has<br>a sound   | This sound is caused by<br>the mechanical operation<br>of the product                                            | This kind of sound is normal,<br>which does not mean that<br>the function of the product fails.                                                |

## Certificate

Product Model : \_\_\_\_\_

Inspector : \_\_\_\_\_

Date of Manufacture : \_\_\_\_\_

Company:

Telephone:

The product is qualified and  
approved to leave the factory

### Warranty Card

#### Precautions for Warranty Conditions

Within one year from the date of purchase, you can enjoy free parts  
replacement and repair services with a valid warranty card or invoice.  
(The replacement of the product shell is not covered by the warranty.  
We promise that the free maintenance includes maintenance labor  
costs and material costs but doesn't includes transportation costs, at  
the same time we do not provide door-to-door maintenance services)  
We are willing to take responsibility for the recycling and reuse of  
products after they are discarded.

Any damages or malfunction caused by improper use by the customer,  
such as an accident caused by supplying wrong voltage or repairing the  
product by personnel not authorized by our company is not covered by  
the warranty.

Any additional services or shipping services that are not covered by the  
standard warranty will be charged accordingly.

Damage caused by improper use, storage and maintenance by  
customers; damage caused by disassembling by maintenance personnel  
not authorized by our company; no warranty card or invoice of our  
company; damage caused by force majeure, our company can provide  
maintenance services, however, the corresponding accessories and  
transportation costs will be charged.

Purchase Model: \_\_\_\_\_

Body Code: \_\_\_\_\_

Date of Purchase: \_\_\_\_\_

Purchase Site: \_\_\_\_\_

Customer Name: \_\_\_\_\_

Contact Number: \_\_\_\_\_

Contact Address: \_\_\_\_\_

Note: due to the continuous upgrading of service, we reserve the  
right to modify the rules without notice.